

Ro-Tap® Sieve Shaker

SV012

Ro-Tap® Sieve Shaker for 200mm & 8 inch diameter sieves

Ro-Tap® Motorised Sieve Shaker for use with 200mm & 8 inch sieves.

The Ro-Tap® shaker accepts up to 6 full-height 200mm or 8 inch sieves or 13 half-height 200mm or 8 inch sieves, plus a full height receiver.

The Ro-Tap Sieve Shaker provides a compact design and aggressive sieving action. This shaker provides 278 oscillations and 150 taps per minute to produce an effective sieving action. It is powered by a 1/4 hp motor and provides a built-in 99 minute step-down timer.



Specification

Sieve sizes accepted	200mm & 8 inch diameter
Timer	99 minute step-down
Motor	1/4 hp
Power Supply	220V 50Hz 1ph
Dimensions	711 x 533 x 635mm
Weight	86kg

Ordering Information

SV012 Ro-Tap® Sieve Shaker for 200mm and 8inch diameter sieves, 220V 60Hz 1ph.